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# SHRINK UNIVERSE

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- 1 Aşiran Serdar, Z., Aktaş Karabay, E., Tallı Parmak, A., & Aksoy, B.(2019). Efficacy of high-intensity focused ultrasound in facial and neck rejuvenation. *Journal of Cosmetic Dermatology*, 19(2), 353-358. DOI: 10.1111/jocd.13008.(Demonstrated significant improvement ) 80% in nasolabial folds, jawline, submental, and neck with appropriate HIFU parameters; recommended moderate energy for safety.(HIFU offers visible improvements) (Efficacy of highintensity focused ultrasound in facial and neck rejuvenation.)
- 2 Hwang, Y., Wan, J., & Yi, K.-H.(2025). Do Different High-Intensity-Focused Ultrasound Frequencies Have Different Effects? A Histological Analysis Correlated With Patients' Subjective Assessments. *Journal of Cosmetic Dermatology*, 24(1), e70069. DOI: 10.1111/jocd.70069.(Histological study using Ultraformer MPT showing increased types I/III collagen and significant fat reduction post-HIFU; clinically 70.6% noted improved facial contour.(Do Different High-Intensity-Focused Ultrasound Frequencies Have Different Effects? A Histological Analysis Correlated With Patients

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(Subjective Assessments.)

- 3 Smoczok, M., Leonik, S., & Bergler-Czop, B.(2022). High-intensity focused ultrasound technology as a non-surgical alternative to face lifting. *Dermatology Review*, 109(2), 130-137. DOI: 10.5114/dr.2022.117984.(Review of HIFU in aesthetics - discusses mechanism of creating thermal coagulation points in SMAS/dermis for tightening and notes fat atrophy as a possible effect and side effect(High-intensity focused ultrasound technology as a non-surgical alternative to face lifting), highlighting the need for proper technique.)
- 4 Park, H., Kim, E., Kim, J., Ro, Y., & Ko, J.(2015). High-intensity focused ultrasound for the treatment of wrinkles and skin laxity in seven different facial areas. *Annals of Dermatology*, 27(6), 688-693. DOI: 10.5021/ad.2015.27.6.688.(Early clinical study in Asian patients showing HIFU improved skin laxity across multiple facial regions; underscores importance of adequate shot delivery, ~420 shots on average, for noticeable tightening.)
- 5 Fritz, K., & Salavastru, C.(2016). Ways of noninvasive facial skin tightening and fat reduction. *Facial Plastic Surgery*, 32(3), 276-282. DOI: 10.1055/s-0036-1584214. (Review of non-invasive modalities including HIFU, explaining how focused ultrasound can target the SMAS for tightening and subcutaneous tissue for fat reduction, enabling facial contour improvements in appropriate patients.)